

Supplementary Table 2: SNP marker used in our assay for eye and hair color prediction with details on PCR primers.

SNP-ID	gene	forward primer	reversed primer	ref.
rs1129038	HERC2	CAGCAGCGACGATTAGATA	ATCACGGCCAGTCAGTCTCT	Eiberg, 2008
rs12913832	HERC2	TGATAGCGTGCAGAACCTTGACA	TAATTCAAAATGCCCCCAAG	Eiberg, 2008
rs1393350	TYR	GGTGAATGATAAACACGAACAGATT	GATGCGTGCATATCCACCAACT	Sulem, 2007
rs7170989	OCA2	TCCCAATTTCGTCTCCTG	TTTCCAGCCAACAAATGAA	Duffy, 2007
rs7495174	OCA2	TGCTGTCGAACCTTCAA	GCTTAGGAAGCAAGGCAAGTT	Duffy, 2007
rs1800407	OCA2	TGAAAGGCTGCCTCTGTTCT	TATGGTCACAGGCGTGAAGA	Duffy, 2007
rs1667394	OCA2	TCCACCATTAAAGACGCAGCAATTCAA	GAGAACTTGAGGTCTCCAAC	Sulem, 2007
rs26722	SLC45A2	TGCCAGCTCTGGATTACG	GCTGGGAAGACTGTTGGGTA	Soejima, 2007
rs16891982	SLC45A2	GGAAAACACGGAGTTGATGC	CCAGAGGTGGAGAACAGAG	Soejima, 2007
rs1426654	SLC24A5	TTCAGCCCTGGATTGTCTC	AATTGCAGATCCAAGGATGG	Soejima, 2007
rs1805005	MC1R	CTGGTGAGCTGGTGGAGA	TCCAGCAGGAGGATGACG	Mengel-From, 2008
rs1805008	MC1R	CTGTCCAGCCTCTGCTTCCT	GGCCACGTGGTCGTAGTAGG	Mengel-From, 2008