

Supplementary Table 2. Point heteroplasmies (PHP) detected in samples used for comparison of libraries produced by 12-cycle and 15-cycle limited-cycle PCR during Nextera® XT library preparation.

Sample ID	PHP	Major/minor allele	Minor allele %	
			12 cycles	15 cycles
MW-004	16519Y	T/C	4.4	4.5
MW-019	16294Y	C/T	3.2	3.7
MW-020	152Y	T/C	5.6	5.7
	9325Y	T/C	5.4	6.2
MW-026	9740Y	C/T	20.9	20.6
	14539R	A/G	29.4	31.0
	16093Y	C/T	14.9	14.1
MW-152	12285Y	T/C	21.6	23.3
MW-236	152Y	T/C	15.6	15.1
	9111Y	T/C	44.6	45.5
	16189Y	T/C	21.7	19.9
MW-290	6899R	G/A	5.3	5.7